# Levenhuk D70L Digital biological microscope

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**levenhuk**   
Zoom&Joy

- 1 Monitor
- 2 LCD screen
- 3 Tubus
- 4 Revolving nosepiece
- 5 Objective lenses
- 6 Arm
- 7 Stage
- 8 Diaphragm disc
- 9 Illumination
- 10 Focus adjustment knob
- 11 Illumination adjustment wheel
- 12 Base
- 13 AC jack

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## General Information

This digital microscope is designed for bright field observations of transparent and solid objects in either transmitted or reflected light. It is well-suited for school studies in biology and may be used for microphotography.

The optics yields sharp and bright views throughout the field of view. The standard kit of this microscope even allows you to edit the microphotography images that you take during research. This microscope is manufactured in accordance with all existing regulations and technical documentation; and was certified to be ready for personal and professional use

## Specifications

|                                   |                                                |
|-----------------------------------|------------------------------------------------|
| Magnification                     | 40x-400x                                       |
| Objective lenses                  | 4x, 10x, 40x                                   |
| Linear field                      | 18 mm                                          |
| Tubus length                      | 160 mm (6.3 in)                                |
| Stage size                        | 95 mm x 95 mm (3.7 x 3.7 in)                   |
| Stage movement range              | 0-15 mm (0-0.6 in)                             |
| Illumination source power supply  | Two AA batteries (not included) or power cable |
| Illumination source (lower/upper) | LED lamp                                       |
| LCD monitor                       | 3.6"                                           |
| Matrix                            | 2 Mpx CMOS (1600x1200, 800x600)                |
| Image format                      | JPG, AVI                                       |
| Interface                         | USB, AV                                        |

## Package

- 4x, 10x, 40x objective lenses
- Diaphragm disc
- LCD monitor with SD card slot
- Microtome
- Forceps
- Brine shrimp hatchery
- Flasks with yeast, brine shrimp eggs, pitch, and sea salt
- 5 ready-to-use specimens and 5 blank slides
- Dust cover
- Plastic case
- User manual and warranty

## Parts of the microscope

**Base:** supports the weight of the microscope and houses the illumination source, electronics and control mechanisms.

**Arm:** holds the base, the stage and the head of the microscope together.

**Rack-and-pinion mechanism:** is mounted on the microscope arm, provides for vertical movements of the stage. Coaxial control knobs are placed on both sides of the arm. Total adjustment range is 15 mm.

**Stage:** is equipped with a diaphragm disc and a specimen holder that can be used to move slides during your observations.

**Locking screw:** to avoid any damage to the objective lens, the microscope arm has a locking screw that acts as an upper limit for stage movements.

**Revolving nosepiece:** houses three objective lenses of the microscope, and allows you to change these lenses smoothly and easily.

**LCD monitor:** instead of an eyepiece this microscope uses a 3.6" LCD monitor.

## Choosing the correct objective lens

Always begin your observations at the lowest magnification possible. This will allow you to center the specimen in your field of view. This alignment of the specimen has to be done carefully, or you might not get the desired view upon changing the magnification.

## Focusing the view

- Turn on the LCD monitor.
- Place a specimen on the stage and use stage clips to hold it in place.
- Turn the revolving nosepiece until 4x objective lens is locked in place. Move the slide manually, to center the desired part of the observed object in the field of view.
- Raise the stage by turning the focus knobs until the cover glass of the slide touches the objective lens. Be careful and try not to damage the slide or the lens.
- Look at the monitor and, if necessary, adjust focus with a fine focus knob until you see a clear image on the screen.

If this procedure is performed properly, you will not have to worry about damaging the optics or the specimen, when switching to a different objective lens. However, additional focusing may still be required.

## Environment conditions

This model performs best under the following conditions: room temperatures from 10°C to 35°C (50°F to 95°F), relative humidity not exceeding 80% at 25°C (77°F). The room should also be void of any dust, acid or base vapors and other active substances

## Battery safety instructions

- Always purchase the correct size and grade of battery most suitable for the intended use.
- Always replace the whole set of batteries at one time; taking care not to mix old and new ones, or batteries of different types.
- Clean the battery contacts and also those of the device prior to battery installation.
- Make sure the batteries are installed correctly with regard to polarity (+ and -).
- Remove batteries from equipment that is not to be used for an extended period of time.
- Remove used batteries promptly.
- Never attempt to recharge primary batteries as this may cause leakage, fire, or explosion.
- Never short-circuit batteries as this may lead to high temperatures, leakage, or explosion.
- Never heat batteries in order to revive them.
- Do not disassemble batteries.
- Remember to switch off devices after use.
- Keep batteries out of the reach of children, to avoid risk of ingestion, suffocation, or poisoning.
- Utilize used batteries as prescribed by your country laws.

## Care and maintenance

- Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS.
- Take necessary precautions when using the device with children or others who have not read or who do not fully understand these instructions.
- After unpacking your microscope and before using it for the first time check for integrity and durability of every component and connection.
- Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center.
- Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws.
- Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special cleaning wipes and special optics cleaning tools from Levenhuk. Do not use any corrosive or acetone-based fluids to clean the optics.
- Abrasive particles, such as sand, should not be wiped off lenses, but instead blown off or brushed away with a soft brush.
- Do not use the device for lengthy periods of time, or leave it unattended in direct sunlight. Keep the device away from water and high humidity.
- Be careful during your observations, always replace the dust cover after you are finished with observations to protect the device from dust and stains.
- If you are not using your microscope for extended periods of time, store the objective lenses and eyepieces separately from the microscope.
- Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fire and other sources of high temperatures.
- When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may heat up during use, and might become a fire hazard.
- Always unplug the microscope from a power source before opening the base or changing the illumination lamp. Regardless of the lamp type (halogen or incandescent), give it some time to cool down before trying to change it, and always change it to a lamp of the same type.
- Always use the power supply with the proper voltage, i.e. indicated in the specifications of your new microscope. Plugging the instrument into a different power outlet may damage the electric circuitry of the microscope, burn out the lamp, or even cause a short circuit.
- Seek medical advice immediately if a small part or a battery is swallowed.

## Troubleshooting

If you are experiencing some problems with your instrument during its service, refer to this section for possible solutions. However, if this doesn't help, you should contact your local service center.

| Problem                                                            | Possible cause                                                                          | Solution                                                                  |
|--------------------------------------------------------------------|-----------------------------------------------------------------------------------------|---------------------------------------------------------------------------|
| Uneven illumination throughout the field of view                   | The objective lens is not locked in place, i.e. optical system is not properly aligned. | Rotate the revolving nosepiece until the lens locks in place.             |
|                                                                    | The diaphragm disc is not locked in place                                               | Rotate the disc until it locks in place                                   |
|                                                                    | One of the lenses is dirty                                                              | Clean the lens from dust or dirt                                          |
| Clean the stage from dust                                          | One of the lenses is dirty                                                              | Clean the lens from dust or dirt                                          |
|                                                                    | The stage is dirty                                                                      | Clean the stage from dust                                                 |
| Reduced quality of views (low resolution, poor contrast)           | One of the lenses is dirty                                                              | Clean the lens from dust or dirt                                          |
| The field of view is clipped                                       | The objective lens is not locked in place                                               | Rotate the revolving nosepiece until the lens locks in place              |
|                                                                    | The diaphragm disc is not locked in place                                               | Rotate the disc to lock the diaphragm in place                            |
|                                                                    | The slide is moving                                                                     | Use stage clips to hold it in place                                       |
| When changing objective lenses, they scrape the slide on the stage | The slide is upside down                                                                | Place it correctly                                                        |
|                                                                    | The cover glass is too thick                                                            | Use a cover glass of standard thickness                                   |
| The illumination system does not work                              | One of the batteries is needs replacing                                                 | One of the batteries is needs replacing Change                            |
| The monitor does not work                                          | The power cable is not plugged in                                                       | Make sure that the power cable is plugged into the socket and the AC jack |

## Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars and other optical products, except for accessories, carry a **lifetime warranty** against defects in materials and workmanship.

Lifetime warranty is a guarantee on the lifetime of the product on the market.

All Levenhuk accessories are warranted to be free of defects in materials and workmanship for **six months** from date of retail purchase.

Levenhuk will repair or replace such product or part thereof which, upon inspection by Levenhuk, is found to be defective in materials or workmanship.

As a condition to the obligation of Levenhuk to repair or replace such product, the product must be returned to Levenhuk together with proof of purchase satisfactory to Levenhuk.

This warranty does not cover consumable parts, such as bulbs (electrical, LED, halogen, energy-saving and other types of lamps), batteries (rechargeable and non-rechargeable), electrical consumables etc.

For further details, please visit our web site: <https://www.levenhuk.com/warranty>

If warranty problems arise, or if you need assistance in using your product, contact the local Levenhuk branch.

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Levenhuk® is a registered trademark of Levenhuk, Inc.
© 2006-2016 Levenhuk, Inc. All rights reserved. 20160428
Levenhuk, Inc. 924-D East 124th Ave. Tampa, FL 33612 USA