### **Europe / Middle-East / Africa**

Microlife AG

Espenstrasse 139

9443 Widnau / Switzerland

Tel. +41 / 71 727 70 30

Fax +41 / 71 727 70 39

Email [admin@microlife.ch](mailto:admin@microlife.ch)

[www.microlife.com](http://www.microlife.com)

### **Asia**

Microlife Corporation.

9F, 431, RuiGang Road, NeiHu

Taipei, 114, Taiwan, R.O.C.

Tel. 886 2 8797-1288

Fax 886 2 8797-1283

Email [service@microlife.com.tw](mailto:service@microlife.com.tw)

[www.microlife.com](http://www.microlife.com)

### **North / Central / South America**

Microlife USA, Inc.

424 Skinner Blvd., Suite C

Dunedin, FL 34698 / USA

Tel. +1 727 451 0484

Fax +1 727 451 0492

Email [msa@microlifeusa.com](mailto:msa@microlifeusa.com)

[www.microlife.com](http://www.microlife.com)

## **microlife BP 3BT0-A**

### **Blood Pressure Monitor for Pregnant Women**

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**microlife®**

# **Blood Pressure Monitor for Pregnant Women**

## Instruction Manual

## **Table of Contents**

### **1. Introduction**

---

- 1.1. Features
- 1.2. Important information about self-measurement

### **2. Important Information on the Subject of Blood Pressure and its Measurement**

---

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?

### **3. Components of your Blood Pressure Monitor**

---

### **4. Using your Blood Pressure Monitor**

---

- 4.1. Inserting the batteries
- 4.2. Using an mains adapter (special accessory)
- 4.3. Cuff connection
- 4.4. Setting the time and date

### **5. Carrying out a Measurement**

---

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Discontinuing a measurement
- 5.6. Memory – recall of the measurements
- 5.7. Memory – cancellation of all measurements

### **6. Error Messages/Malfunctions**

---

### **7. Care and Maintenance, Recalibration**

---

### **8. Guarantee**

---

### **9. Technical Specifications**

---

### **10. [www.microlife.com](http://www.microlife.com)**

---

## 1. INTRODUCTION

---

### 1.1. Features

Your blood pressure monitor is a fully automatic, digital blood pressure measuring device for use on the upper arm. This monitor enables very fast and reliable measurement of the systolic and diastolic blood pressure as well your pulse by use of the oscillometric method.

Up to 20 % of women develop hypertension during pregnancy. This disease (pre-eclampsia or «toxemia») affects the final outcome of about 5 % of all pregnancies. Pre-eclampsia is a disorder that occurs only during pregnancy. It can be recognized by a clear increase in blood pressure and high protein levels in the urine. The measurement precision of this device has been proven clinically in a comprehensive study carried out by Professor Andrew Shennan at St. Thomas' Hospital, London. This study proved the reading reliability of this monitor during pregnancy and in the presence of pre-eclampsia, making it the first such home-use monitor that is clinically validated for such use and especially suitable for pregnant women. Of course, this device can be used also after the pregnancy or by other family members.

### Attention!

#### 1.2. Important information about self-measurement

- Self-measurement means **Control**, not diagnosis or treatment. Your values must always be discussed with your doctor. **Never alter the dosages of any medication without direction from your doctor.**
- The pulse reading is **not** suitable for checking heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

#### Electromagnetic interference:

The device contains sensitive electronic components. Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

## 2. Important Information on the Subject of Blood Pressure and its Measurement

---

### 2.1. How does high/low blood pressure arise?

Your blood pressure level is determined in the circulatory center of your brain. Your nervous system allows your body to adapt or alter blood pressure in response to different situations. Your body alters your pulse and the width of blood vessels through changes in muscles in the walls of blood vessels.

Your blood pressure reading is highest when the heart pumps or ejects blood – this is called your Systolic Blood Pressure. Your blood pressure reading is lowest when the heart rests (in-between beats) – this is called your Diastolic Blood Pressure.

**Note:** Maintain blood pressure values within a «normal» range(s) in order to prevent particular diseases.

### 2.2. Which values are normal?

Blood pressure is too high if at rest, the diastolic pressure is above 90 mmHg and/or the systolic blood pressure is over 140 mmHg. If you obtain readings in this range, consult your doctor immediately. High blood pressure values over time damage blood vessels, vital organs such as the kidney and even your heart.

When blood pressure values are too low, i.e. systolic values under 100 mmHg and/or diastolic values under 60 mmHg, consult your doctor.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to independently alter drug doses prescribed by your doctor.**

Table for classifying blood pressure values (units mmHg) according to World Health Organization:

| Range                 | Systolic Blood Pressure | Diastolic Blood Pressure | Measures                        |
|-----------------------|-------------------------|--------------------------|---------------------------------|
| Hypotension           | lower than 100          | lower than 60            | Consult your doctor             |
| Optimal range         | between 100 and 120     | between 60 and 80        | Self-check                      |
| Normal range          | between 120 and 130     | between 80 and 85        | Self-check                      |
| High-normal range     | between 130 and 140     | between 85 and 90        | Consult your doctor             |
| Mild hypertension     | between 140 and 160     | between 90 and 100       | Consult your doctor             |
| Moderate hypertension | between 160 and 180     | between 100 and 110      | Consult your doctor             |
| Severe hypertension   | higher than 180         | higher than 110          | Consult your doctor immediately |

**Attention!** Adjusted values **during pregnancy** (units mmHg):

| Range               | Systolic Blood Pressure | Diastolic Blood Pressure | Measures                        |
|---------------------|-------------------------|--------------------------|---------------------------------|
| Normal range        | lower than 140          | lower than 90            | Self-check                      |
| Hypertension        | higher than 140         | higher than 90           | Consult your doctor             |
| Severe hypertension | higher than 160         | higher than 100          | Consult your doctor immediately |

*Reference: Prof. A.H. Shennan, St. Thomas Hospital, London*

### 2.3. What can be done, if regular increased/low values are obtained?

- Please consult your doctor.
- Increased blood pressure values (various forms of hypertension) over time are associated with considerable risks to health. Blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). This can result in a deficient supply of blood to important organs (heart, brain, muscles). Additionally, the heart will become structurally damaged.
- There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. Secondary hypertension can cause organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- There are lifestyle changes you can make to prevent and reduce high blood pressure. These measures must be part of a healthy lifestyle and include:

### A) Eating habits

- Strive for a normal weight as prescribed by your doctor. Reduce overweight!
- Avoid excessive consumption of common salt. (Many «packaged foods» contain high levels of salt).
- Avoid fatty foods. (Packaged foods are frequently high in fats).

### B) Previous illnesses

Consistently follow any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus or sugar diabetes)
- Fat metabolism disorder
- Gout

### C) Habits

- Eliminate smoking
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee, tea, chocolate, etc.)

### D) Physical conditioning


- After a preliminary medical examination, exercise regularly.
- Choose sports which require endurance and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise activities. Your physician will help you develop an exercise routine that is appropriate for you.

## 3. Components of the Blood Pressure Monitor

---


The illustration shows the blood pressure monitor, consisting of:

### a) Measuring unit

### b) Upper arm cuff

M-size-cuff (22 - 32 cm) and  
L-size-cuff (32 - 42 cm)


## 4. Using your Blood Pressure Monitor

### 4.1. Inserting the batteries


Insert batteries immediately after unpacking the device. The battery compartment is located on the bottom of the device (see illustration).

- Remove cover as illustrated
- Insert the batteries (4 x size AA 1.5V), thereby observing the indicated polarity.
- If the battery warning appears in the display, the batteries are empty and must be replaced.

### Attention!


- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AA» Long-Life or Alkaline 1.5 V Batteries. The use of 1.2 V Accumulators is not recommended.
- If the blood pressure monitor is left unused for long periods, please remove the batteries from the device.

### 4.2. Using a mains adapter (special accessory)

This blood pressure monitor can be operated with the Microlife mains adapter (output 6 V DC / 600 mA, DIN plug).

- Plug the plug into the socket at the back of the instrument.
- Plug the mains adapter into a 230 V or 110 V power socket. Test that power is available by pressing the O/I button.


### Note:

- No power is taken from the batteries while the mains adapter is connected to the instrument.
- If the mains voltage is interrupted during the measurement (e.g. by accidental removal of the mains adapter from the mains socket), the instrument must be reset by removing the plug from its socket and re-inserting the connections.
- Please consult a specialist dealer if you have questions relating to the mains adaptor.

### 4.3. Cuff connection

Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.



#### 4.4. Setting the time and date

This blood pressure monitor automatically records the time and date of each measurement. This is very important information as blood pressure normally varies over the course of a day.

1. After new batteries have been inserted the year is blinking in the display. You can adjust the year by pressing the «M» button. Press the Time button to confirm and switch to month setting.

2. The correct month can be entered by the «M» button. Press the Time button to confirm and switch to day setting.

3. Please follow the instructions above to set the day, hour and minutes.
4. Once the minutes have been set and the Time button is pressed, the date and time will be displayed on the screen.
5. If you want to change the date and time again, please hold the Time button down for about 3 seconds until the year number starts to flash. Now you can enter the new values as described above.

### 5. Carrying out a Measurement

---

#### 5.1. Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try to find time to relax by sitting in an armchair in a quiet atmosphere for about 5 minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Take measurements regularly at the same time of day, since blood pressure changes during the course of the day.

#### 5.2. Common sources of error

**Note:** Comparable blood pressure measurements always require the same conditions!

These are normally always quiet conditions.




- All efforts by the patient to support their arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (higher) than the heart, an false reading will be obtained. Each 15 cm difference in height results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is

of extraordinary importance. The cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please contact your dealer. **Note: Only use clinically approved Microlife Cuffs!**

- A loose cuff or a sideways protruding air pocket causes false measurement values.

### 5.3. Fitting the cuff

- Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and that the rubber tube leaves the cuff on the inner side of the arm.
- Tighten the free end of the cuff and close the cuff with the Velcro.
- The cuff should be snug on your upper arm, but not too tight. Any clothing which restricts the arm (e.g. pullover) must be taken off.
- Lay your arm on a table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked.


#### Note:

If it is not possible to fit the cuff to the left arm, it can also be placed on the right. However all measurements should be made using the same arm.

### 5.4. Measuring procedure

After the cuff has been appropriately positioned, the measurement can begin:


- Press the START button. The pump begins to inflate the cuff. The rising pressure in the cuff is shown in the display.
- After reaching the inflation pressure, the pump stops and the pressure gradually falls. The cuff pressures are displayed. In case that the inflation pressure is not sufficient, the monitor automatically re-inflates to a higher level.
- When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep is heard for every heartbeat.
- A longer beep is sounded when the measurement has been completed. The systolic, and diastolic blood pressures and pulse rate now appear in the display.
- The measurement results are displayed, until you switch the device off. If no button is pressed for 5 minutes, the device switches automatically off, to save the batteries.

### 5.5. Discontinuing a measurement


If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the START button can be pressed at any time.

The device immediately lowers the cuff pressure automatically.

### 5.6. Memory - recall of measurements

This blood-pressure monitor automatically stores the last 30 measurement values. By pressing the «M» button, the last Measurement, as well as previous measurements can be displayed one after the other.

(MR30: Value of the last measurement)

(MR29: Value of the measurement before MR 30)



### 5.7. Memory- cancellation of all measurements

#### Attention!


Before you delete all readings stored in the memory, make sure you will not need refer to the readings at a later date.

In order to delete all stored readings, depress the «M» button for at least 7 seconds, the display will show the symbol «CL» and 3 short beep sounds will be heard to indicate deletion of stored readings.

## 6. Error Messages/Malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.  
(Example: Error no. 1)

| Error No. | Possible cause(s)                                                                                                                                                                                                                   |
|-----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ERR 1     | The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. The tube may have become unplugged after the systolic blood pressure was measured.<br>Further possible cause: No pulse has been detected. |
| ERR 2     | Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artefact).                                                                                                          |

| Error No. | Possible cause(s)                                                                                                                                                                                                      |
|-----------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ERR 3     | Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed.                                                                                                           |
| ERR 5     | The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings. |
| HI        | The measured pulse is higher than 200 beats/min.                                                                                                                                                                       |
| LO        | The measured pulse is lower than 40 beats/min.                                                                                                                                                                         |

**Other possible malfunctions and possible solutions** – If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

| Malfunction                                                                                                      | Remedy                                                                                                                                                                                                                                                                                                                                                                                                       |
|------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The display remains empty when the instrument is switched on although the batteries are in place.                | <ol style="list-style-type: none"> <li>1. Check whether the batteries are installed with the correct polarity and correct if necessary.</li> <li>2. If the display is unusual, remove the batteries and then exchange them for new.</li> </ol>                                                                                                                                                               |
| The pressure does not rise although the pump is running.                                                         | <ul style="list-style-type: none"> <li>• Check the connection of the cuff tube and connect properly if necessary.</li> </ul>                                                                                                                                                                                                                                                                                 |
| The device frequently fails to measure the blood pressure values, or the values measured are too low (too high). | <ol style="list-style-type: none"> <li>1. Fit the cuff correctly on the arm.</li> <li>2. Before starting the measurement make sure that the cuff is not fitted too tightly and that there is no rolled-up sleeve exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary.</li> <li>3. Measure the blood pressure again in complete peace and quiet.</li> </ol> |
| Every measurement results in different values although the device functions normally and the values displayed.   | <ul style="list-style-type: none"> <li>• Please read the following information and the points listed under «Common sources of error». Repeat the measurement.</li> </ul>                                                                                                                                                                                                                                     |
| The blood pressure values measured differ from those measured by the doctor.                                     | <ul style="list-style-type: none"> <li>• Record the daily development of the values and consult your doctor about them.</li> </ul>                                                                                                                                                                                                                                                                           |

### Further information

The level of blood pressure is subject to fluctuations even in healthy people. It is important to compare measurements taken under the same conditions and at the same time of day. (Quiet conditions)!



If you have any questions regarding the use of this blood pressure monitor, please ask your dealer or pharmacist for the Microlife Service representative in your country. The Microlife Service Team will be happy to help you. **Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all guarantee claims!

## 7. Care and Maintenance, Recalibration

---

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle the cuff carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Carefully remove spots on the cuff with a damp cloth and soapsuds. **Do not wash the cuff in a clothes or dishwasher!**
- d) Handle the tube carefully. Avoid stress from pulling, pinching and sharp edges.
- e) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the device!** This can negatively impact the calibration! (accuracy)
- g) The batteries should be removed, if the device is not going to be used for an extended period of time. Do not continue to use batteries beyond the expiry date indicated.

### Periodical re-calibration

Sensitive measuring devices must be checked for accuracy from time to time. We recommend a periodical inspection of your device by an authorized microlife dealer **every 2 years**. Your factory authorized Microlife dealer would be pleased to provide more information regarding calibration.

## 8. Guarantee

---

This blood pressure monitor is **guaranteed for 3 years** from date of purchase. This guarantee includes the instrument and the cuff. The guarantee does not apply to damage caused by improper handling, damage from leaking batteries, accidents, not following the operating instructions or alterations made to the instrument by third parties.

The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

Name and company address of the responsible dealer:

## 9. Technical Specifications

---

|                                     |                                                                                     |
|-------------------------------------|-------------------------------------------------------------------------------------|
| <b>Weight:</b>                      | 460 g (with batteries)                                                              |
| <b>Size:</b>                        | 115 (W) x 182 (L) x 76 (H) mm                                                       |
| <b>Storage temperature:</b>         | -20 to +50°C                                                                        |
| <b>Humidity:</b>                    | 15 to 90% relative humidity maximum                                                 |
| <b>Operation temperature:</b>       | 10 to 40°C                                                                          |
| <b>Display:</b>                     | LCD-Display (Liquid Crystal Display)                                                |
| <b>Measuring method:</b>            | oscillometric                                                                       |
| <b>Pressure sensor:</b>             | capacitive                                                                          |
| <b>Measuring range:</b>             |                                                                                     |
| <b>SYS/DIA:</b>                     | 30 to 280 mmHg                                                                      |
| <b>Pulse:</b>                       | 40 to 200 beats per minute                                                          |
| <b>Cuff pressure display range:</b> | 0–299 mmHg                                                                          |
| <b>Memory:</b>                      | Automatically stores the last 30 measurements                                       |
| <b>Measuring resolution:</b>        | 1 mmHg                                                                              |
| <b>Accuracy:</b>                    | Pressure within $\pm 3$ mmHg                                                        |
| <b>Pulse</b>                        | $\pm 5$ % of the reading                                                            |
| <b>Power source:</b>                | a) 4 x 1.5V batteries; LR6, size AA<br>b) mains adapter 6V DC 600mA (optional)      |
| <b>Cuffs:</b>                       | M-size-cuff (22 - 32 cm) and<br>L-size-cuff (32 - 42 cm)                            |
| <b>Reference to standards:</b>      | EU directive 93/42/EEC<br>NIBP - requirements: EN 1060-1 /-3 /-4 & ANSI / AAMI SP10 |

Technical alterations reserved!

## 10. [www.microlife.com](http://www.microlife.com)

---

Detailed user information about our products as well as services can be found at [www.microlife.com](http://www.microlife.com).

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